


DISINFECTANT
PT2 & PT4 BIOCID




UMONIUM³⁸ MASTER SPRAY

✓ FORMULA FREE FROM

- CMR* COMPONENT
- ENDOCRINE DISRUPTOR

FLASH DISINFECTION
AFTER 30 SEC.



“ Disinfects within 30 seconds
of contact time ”

Ultra-fast-acting disinfectant covering the whole spectrum, **UMONIUM³⁸ MASTER SPRAY** kills all germs present on treated surfaces with a single spray. Its performance is particularly guaranteed in situations where many disinfectants are powerless. It easily eliminates organic residues such as traces of blood, urine and faeces. Simple and ergonomic, the double nozzle function projects an ideal foam for cleaning or sprays out a fine jet for optimal disinfection.

UMONIUM³⁸ MASTER SPRAY is a biocidal product. Use biocides safely. Read the label and product information before use.

* Carcinogenic, Mutagenic or Reprotoxic

 **laboratoire
HUCKERT'S**
INTERNATIONAL
mastered disinfection





UMONIUM³⁸ MASTER SPRAY

1L amounts to 1000 sprays and treats an area of 200 m².



✓ FORMULA FREE FROM

- CMR* COMPONENT
- ENDOCRINE DISRUPTOR

FIELDS OF APPLICATION

Low and high surfaces, materials, equipment, furniture and sanitary facilities, medical transport vehicles, or transportation of the deceased, hospital waste, ventilation systems (filters and ducts), healthcare premises, hospital laundry, etc.

COMPOSITION

Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-16)). 3.3g/L.

REGULATORY STANDARDS

- Quality and traceability guaranteed by ISO 9001:2015
- Complies with biocide regulation 528/2012

INSTRUCTIONS FOR USE

Cleaning/Pre-disinfection: Spray in the "foam" position. Brush if necessary. Rinse.

Disinfection: Use (in fine position) on clean surfaces. Ensure a minimum contact time of 15 minutes. Rinse and wipe the surface. For surfaces in contact with food products, rinse twice.

PROPERTIES


- PT2 & PT4 biocide
- Detergent- disinfectant with a broad spectrum of activity
- Non-toxic: LD₅₀ (rat 5840 mg/kg), no dangerous vapours
- Non-corrosive: pH-neutral
- Ready-to-use product, non-irritating for the skin (no transdermal toxicity)
- Free from aldehyde, peroxide and biguanide, and other carcinogenic, mutagenic, or reprotoxic components
- Non-coagulant
- Non-flammable
- Bactericide (including MRSA and other antibiotic resistant strains), fungicide, mycobactericide and virucide (including human coronavirus)
- 100% compatible on fibre optic, rubber, polycarbonate, acrylic compound, glass, Pyrex, PU, PVC, HDPE, PET, neoprene, latex, silicone, paints and stainless steel 410
- Effective, fast-acting and efficient (dissolves blood)
- Effective in the presence of organic materials
- Acts within 30 seconds
- No residue after rinsing
- Stability: 36 months date of manufacture

PROCEDURE FOR USE


STEP 1: CLEANING/PRE-DISINFECTION



Set the position to "foam"



Spray



Use a clean cloth




Wash carefully



Rinse


STEP 2: DISINFECTION



Set the position to fine spray



Spray



Leave to act for a minimum of 15 min.



Rinse*



Wipe with a clean cloth

* For surfaces in contact with food products, rinse twice.

TECHNICAL DETAILS

| Standards | Target organisms | Contact time (min.) | |
|-----------|---|---------------------|-----------------|
| | | Condition clean | Dirty condition |
| EN 1276 | <i>S. aureus</i> , <i>E. hirae</i> , <i>E. coli</i> , <i>P. aeruginosa</i> | 5 | 10 |
| EN 13727 | <i>S. aureus</i> , <i>E. hirae</i> , <i>P. aeruginosa</i> | 30 sec. | 1 |
| | <i>E. faecalis</i> (vancomycin-resistant), <i>E. coli</i> (carbapenem-resistant), <i>K. pneumoniae</i> (carbapenem-resistant), <i>S. aureus</i> (methicillin-resistant) | 30 sec. | - |
| | <i>S. aureus</i> , <i>E. hirae</i> , <i>E. coli</i> , <i>P. aeruginosa</i> | 10 | 15 |
| | <i>S. typhimurium</i> , <i>E. cloacae</i> , <i>L. brevis</i> | - | 15 |
| EN 16615 | <i>S. aureus</i> , <i>E. hirae</i> , <i>E. P. aeruginosa</i> | 5 | - |
| EN 1650 | <i>C. albicans</i> | 10 | 10 |
| EN 13624 | <i>C. albicans</i> | 5 | 10 |
| EN 13697 | <i>C. albicans</i> | 10 | 15 |
| EN 16615 | <i>C. albicans</i> | 1 | - |
| EN 14476 | Norovirus, Adenovirus, Poliovirus | 10 | 10 |
| | Human coronavirus, Vaccinia virus | 5 | - |
| EN 14348 | <i>M. avium</i> , <i>M. terrae</i> | 10 | 10 |
| EN 14563 | <i>M. avium</i> , <i>M. terrae</i> | 10 | 10 |



| ITEM No. | NAME | BOX |
|----------|--|------------|
| PF 10134 | 250mL spray bottle | 30 x 250mL |
| PF 10130 | 1L spray bottle | 12 x 1L |
| EMB 5206 | 1L Stainless steel bottle bracket | 1 item |
| EMB 0099 | 100 towels UMONIUM ³⁸ Wipes | 12 items |

* Carcinogenic, Mutagenic or Reprotoxic

FP00182 EN - Ed. 3.5 - 12/2023 - Product photos are non-binding.



UMONIUM³⁸ is a registered trademark of Laboratoire Huckert's International