

DERMATOLOGICAL  
SOAP



# PHYTOGEL®

By HUCKERT'S

nature

INGREDIENTS  
95.1%  
OF NATURAL ORIGIN\*



“ A formula that respects  
your skin ”



## DERMATOLOGICAL SOAP FOR DRY TO VERY DRY HANDS


Frequent washing with soap can irritate and damage your skin (flaking, redness, skin spots, etc.). A true soap product, PHYTOGEL By HUCKERT'S® **nature** is also enriched with eucalyptus, thyme and chamomile essential oils. Naturally scented, it respects the skin barrier and protects the driest skin. PHYTOGEL By HUCKERT'S® **nature** is your daily partner in the fight against skin damage.



\* According to ISO 16128

laboratoire  
**HUCKERT'S**  
INTERNATIONAL  
*mastered disinfection*





1 x 400-mL bottle =  
167 doses

### FIELDS OF APPLICATION

- Pre-operative skin preparation
- Frequent use
- Washing dry to very dry skin
- Bio-decontaminating
- Particularly recommended to limit the appearance of chapping

### TECHNICAL DETAILS

- Thick gel, a translucent golden colour
- Characteristic plant odour
- Highly concentrated formula
- Natural perfume (essential oils)
- Dermatologically tested

### COMPOSITION

AQUA, GLYCERIN, OLEIC ACID, AMINOMETHYL PROPANOL, COCONUT ACID, PPG-5-LAURETH-5, COCOAMIDOPROPYL BETAINE, EUCALYPTUS GLOBULUS LEAF OIL, THYME FLOWER/LEAF OIL, HYDROLYSED SOY PROTEIN, PEG-8 CAPRYLIC/CAPRIC GLYCERIDES, CAPRYLYL/CAPRYC ACID, SODIUM CHLORIDE, CHAMOMILLA RECUTITA (MATRICARIA) FLOWER OIL, 2-AMINOBTANOL, SODIUM BENZOATE, BENZOIC ACID, SORBIC ACID, LIMONENE\*\*, LINALOOL\*\*

### REGULATORY STANDARDS

- Quality and traceability guaranteed by ISO 9001:2015 and EN ISO 13485:2016

## PROCEDURE FOR USE

### HAND WASHING



Apply to dry skin



Massage for 15 to 30 sec.



Rinse

### CARE FOR DRY HANDS



Apply to dry skin



Massage for 15 to 30 sec.



Leave on for 1 to 2 min.



Rinse



ITEM No.	NAME	BOX
PF 20112	5-L container with pump	4 x 5 L
EMB 5101	Stainless steel stand for 5L	1 item
PF 20809	400-mL bottle	15 x 400 mL
EMB 6100	Stainless steel 400-mL bracket	1 item
PF 20811	100-mL pocket-sized bottle	28 x 100 mL
EMB 1253	Automatic wall dispenser	1 item
EMB 1252	Manual wall dispenser	1 item
PF 20109	Refill 1-L manual and automatic dispenser	11 x 1L

\* According to ISO 16128

\*\* Naturally present in plant extracts

FP00072 EN - Ed. 3.4 - 01/2023 - Product photos are non-binding.

